



healthkon<sup>®</sup>  
anytime · anywhere · healthcare

# Aarohan

Training Guide-  
System Requirements

## Table of Contents

1. Prerequisites .....	2
2. System Requirements .....	2
3. Browser Settings .....	3

### 1. Prerequisites

This section describes the recommended System requirements and browser settings for successful transactions on the Healthkon Aarohan Platform.

### 2. System Requirements



- Laptop with minimum 2 GB of RAM
- Laptop with Integrated speakers and video camera
- External Mouse
- External Video camera: minimum 5 Megapixels
- Audio port to connect headphones / external speakers
- Integrated Bluetooth
- Windows 2000 / 7 or 10
- Latest Chrome / Firefox browser
- Built-in Wi-Fi / LAN Port
- 3 USB Ports (Preferably USB 3.0)
- HDMI Port to connect to TV/Monitor if required
- LAN / USB Laser Printer with a duty cycle of 15000/Month  
(400 Patients \* 20 working days \* 2 Pages) - Recommended HP or Epson
- Internet connectivity through landline / Cable / Other with Wi-Fi Router
- UPS to take on the device load during power outages
- Printer set up and configured on the Node Admin PC for printing  
consultation report or other documents as required






### 3. Browser Settings

When using the Video Consultation window, allow the use of, Location, Microphone and Camera.



1. On the browser, click **Settings**, scroll down to **Privacy and security > Site Setting > Permissions > Camera**. Verify that the website URL, <https://aaroan.healthkon.com> is allowed access to the Camera.



2. On the browser, click **Settings**, scroll down to **Privacy and security > Site Setting > Permissions > Camera**. Verify that the website URL, <https://aaroan.healthkon.com> is allowed access to the Microphone.

